



IQ MEDICAL SPECIALIST PRODUCT SPECIFICATION

	Dual Bariatric	Dual Summit	Dual Summit Bariatric	Theraflo (LOW AIR LOSS)	Theraflo Bariatric (LOW AIR LOSS)	Dual Bariatric Cushion
Product Order Code	SHCAR-2020	SHCAR-2023	SHCAR-2009	SHCAR-0147	SHCAR-2466	SHNHS-2265
Risk Category	Very High Cat ≤ 4	Very High Cat ≤ 4	Very High, Cat ≤ 4	Very High, Cat ≤ 4	Very High, Cat ≤ 4	Very High, Cat ≤ 4
Max Patient Weight	200kg	300kg	300kg	450kg	450kg	350kg
Cells	20	21	21	20	20	7
Cell Cycle	1 in 2 (10 Minutes)	1 in 2 (10 Minutes)	1 in 2 (10 Minutes)	1 in 2 (10 Minutes)	1 in 2 (10 Minutes)	1 in 2 (10 Minutes)
Foam Underlay	✓	x	x	✓	✓	x
Width (mm)	1200	900	1200	900	1210	750
Length (mm)	2000	2000	2000	2000	2000	500
Depth (mm)	250	250	250	250	250	76
Weight Setting	✓	x	x	✓	✓	x
Comfort Control	x	✓	✓	x	x	✓
Multiple Mode (1)	✓	✓	✓	✓	✓	x
Single Mode (2)	x	x	x	x	x	✓
Cable Management	✓	✓	✓	x	x	x
Pump Lockout	✓	✓	✓	✓	✓	x
Power Failure Alarm	✓	✓	✓	✓	✓	x
Pump Failure Alarm	✓	✓	✓	✓	✓	x
Low Pressure Alarm	✓	✓	✓	✓	✓	✓
Ease of Transport Mode	✓	✓	✓	✓	✓	✓
CPR Emergency Deflation	✓	✓	✓	✓	✓	N/A
Vapour-permeable, high-comfort, Dartex® cover	✓	✓	✓	✓	✓	x
Standard 'breathable' comfort cover	x	x	x	x	x	✓
Weight of Product	21kg	20kg	23kg	19kg	23kg	4.5kg
3 YEAR WARRANTY (3)	Yes	Yes	Yes	Yes	Yes	No

(1) Multiple Mode = Dynamic, Constant Low Pressure, and Maximum Inflate
 (2) Single Mode = Dynamic mode only
 (3) Applicable to purchases made after 1st June 2017
 Please Note all specifications apply to 2017 Product Models