

iQ MEDICAL HYBRID SPECIFICATION	Core Hybrid FoC System	Stellar Hybrid FiC System
Product Order Code: Mattress	IQM-0101	IQM-0201
Product Order Code: Pump	IQM-0100	IQM-0200
Risk Category	Very High Cat ≤4	Very High Cat ≤4
Max Patient Weight: Static	159kg	159kg
Max Patient Weight: Powered	250kg	250kg
Cells	10	14
Cell Cycle	1 in 2 (10 minutes)	1 in 2 (10 minutes)
Foam Over Cell (FoC)	✓	✗
Foam Inside Cell (FiC)	✗	✓
Width (mm)	900	900
Length (mm)	2000	2000
Depth (mm)	150	150
Weight of Mattress	14.0kg	14.0kg
Weight of Pump	1.7kg	1.7kg
24hr care with 'Cushion' connector	✓	✓
Power Failure Alarm	✓	✓
Removable Cover	✗	✓
Easy Grip Handles	✗	✓
Cable Management	✗	✓
CPR Emergency Deflation (Pump)	✓	✓
Vapour-permeable, high-comfort, Dartex® cover	✓	✓
Profiling 'flexi-foam' edges	✓	✓
Specifically designed foam sections for head, body and feet	✓	✓
PUMP WARRANTY	2 years	2 years
COVER WARRANTY	4 years	4 years
FOAM WARRANTY	8 years	8 years
INSERTS WARRANTY	2 years	2 years

iQ Medical

Unit 1a | The Wharf | Stretton Under Fosse | Rugby | CV23 0PR

Telephone: 01788 833 989 | Fax: 0845 127 50 50

Email: sales@iQMedical.co.uk | Web: iQMedical.co.uk